
NGC System

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The background features a light blue gradient on the left side. On the right, there are several overlapping geometric shapes: a vertical orange bar, a vertical blue bar, and a large green shape with rounded corners. A horizontal blue bar crosses the vertical bars. A small brown square is located at the intersection of the horizontal blue bar and the vertical orange bar.

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

BioFrac: especificaciones

The background features a white starburst shape in the upper left quadrant. A thick blue horizontal bar spans the width of the slide, intersecting a vertical orange bar. To the right, there are vertical bars in blue and green, with a dark green rounded rectangular shape overlapping the blue bar and the green bar.

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New Method

Evaluation

Open a run, analyzed run or trace comparison

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
Run Log

Time	Description
15:15:05	Zero UV Baseline 280 nm
15:15:31	Set Sample Inject Valve position - System Pump Inject Loop
15:17:19	Set Sample Inject Valve position - Manual Load Loop / System Pump to Column
15:18:01	Start Fraction Collection - Fraction Size 1 ml

Run Log Run Queue

Pantalla de Control del Sistema

System Pump F10

Flow Control

Flow Rate [0.001-10] 1.000 ml/min

Mode Isocratic Gradient

Start %B 0 End 50

Duration 1 min

Pressure Limits

Lower 0.0 MPa Upper 25.2 MPa

Control Flow to avoid overpressure

Stop Hold

Sample Inject Valve

- Manual Load Loop System Pump to Column
- System Pump Inject Loop
- System Pump Waste

3.1 mS/cm
25.0 °C

UV Con

Time	Description
15:15:05	Zero UV Baseline 280
15:15:31	Set Sample Inject Valve position - System Pump Inject Loop
15:17:19	Set Sample Inject Valve position - Manual Load Loop / System Pump to Column
15:18:01	Start Fraction Collection - Fraction Size 1 ml

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BIO-RAD

Home System Control **Method Editor**

New Open Open Template Save Save As Save As Template Start Run

Template Name: Anion Exchange

Phase Library | **Method Outline** | **Phase Parameters**

Method Settings

Equilibration

Sample Application

Column Wash

Elution

System CIP-Storage

Column Preparation-Activation

Column Performance Test

System Preparation

Create New Phase

Standard Custom

Fluidic Scheme: NGC Quest

Column Selection: Anion Exchange

Column Type: ENrich Q, 5/50 mm

Column Volume: 0.99

Max Pre-C: Custom

Single Wave: 255 nm, 280 nm

Fraction Collection: Rack Type: F1 (12-13 mm x 100 mm tubes)

Flow: Flow Rate: 1.000 [0.001-10] ml/min

Buffer Selection: Manually Prepared Buffer

Inlet A: Buffer A

Inlet B: Buffer B

Run Name

Notes

Configure

Control the flow to avoid overpressure



Detect end of buffer with Air Sensor

Pump Head Type: F10, F100

Name Buffers

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- ... GSTPrep FF 16/10
- ... GSTrap 4B, 1 ml
- ... GSTrap 4B, 5 ml
- ... GSTrap FF, 1 ml
- ... GSTrap FF, 5 ml
- ... GSTrap HP, 1 ml


Durante la corrida...

File View Manual Tools Help Debug

Home System Control Method Editor

BIO-RAD



Method Run - Run AffTest 01 (01/03)

Stop Run Pause Run Hold Step

Fluidic Scheme: AKTA Kryptonite

Affinity Native IMAC Test (00:20:39)

Method Control ? Flow Rate ? UV Monitor ? Fraction Collector ?

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File View Analysis Help

Home System Control Evaluation

Analysis of SP-JH cation... x Analysis of SP-JH cation... x

Add Run... Remove Run... Save Analysis Stack Overlay Peak Integration

SP-JH cation 6/13/2012 12:12:23 PM

A (mAU)

Volume (ml)

Peak Number	Start (ml)	End (ml)	Retention Time (min)	Retention Volume (ml)	Height (mAU)	Area (ml*mAU)	Relative Area
1	-0.92	-0.65	-0.84	-0.83	12.24	1.30	
2	0.52	1.37	0.60	0.60	362.29	67.06	1
3	6.85	7.50	7.32	7.32	34.36	13.12	
4	7.50	8.42	7.88	7.88	464.00	110.56	1
5	8.47	9.10	8.79	8.80	23.80	6.80	
6	9.50	10.85	10.13	10.13	679.54	231.99	3
7	12.62	13.23	12.99	13.00	969.72	135.95	2
8	13.23	13.73	13.33	13.33	74.34	13.43	

Runs/Traces Peaks

Run: SP-JH cation 6/13/2012 12:12:23 PM, Type: λ

Peak Integration

Auto Integration Manual Integration

Traces

SP-JH cation 6/13/2012 12:12:23 PM

λ Automatic

Baseline Parameters

By Best Fit 8.00 (1-10)

By Offset

Peak Parameters

Slope: 10.00 (1-100)

Sensitivity: Medium

Peak Filtering

Min Height: 1.0 mAU

Min Width: 0.001 ml

Size: 20 N-largest

Range

Start: -0.96 ml

End: 24.53 ml

Integrate Delete

Peak Integration

Auto Integration Manual Integration

Traces

SP-JH cation 6/13/2012 12:12:23 PM

λ Automatic

Baseline Parameters

Peak Parameters

Peak Filtering

Default Parameters: Reset

Integrate Delete

User: pjohnso1

Simulating Instrument: NGC

Análisis comparativo

File View Analysis Help


BIO-RAD

Home System Control Evaluation

Trace Comparison 1 * x

Add Run... Remove Run... Save Analysis Stack Overlay Peak Integration

Overlay Trace Comparison 1

Runs/Traces

Show	Type	Color	Min Y	Max Y	Units	Method	Date	Technique
<input checked="" type="checkbox"/>	Run: SP-JH cation 6/13/2012 11:39:36 AM							
<input checked="" type="checkbox"/>	λ	Red	0	2000	mAU	Cation Exchange	6/13/2012 11:39 AM	Undefined
<input type="checkbox"/>	Conductivity	Blue	0	500	mS/cm	Cation Exchange	6/13/2012 11:39 AM	Undefined
<input type="checkbox"/>	%B	Black	0	110	%	Cation Exchange	6/13/2012 11:39 AM	Undefined
<input type="checkbox"/>	System Pressure	Blue	0	500	psi	Cation Exchange	6/13/2012 11:39 AM	Undefined
<input type="checkbox"/>	PreCol Pressure	Green	0	500	psi	Cation Exchange	6/13/2012 11:39 AM	Undefined
<input type="checkbox"/>	ΔCol Pressure	Red	0	500	psi	Cation Exchange	6/13/2012 11:39 AM	Undefined
<input type="checkbox"/>	Temperature	Black	0	100	°C	Cation Exchange	6/13/2012 11:39 AM	Undefined
<input type="checkbox"/>	Flow Rate	Purple	0	100	ml/min	Cation Exchange	6/13/2012 11:39 AM	Undefined
<input checked="" type="checkbox"/>	Run: SP-JH cation 6/13/2012 12:12:23 PM							
<input checked="" type="checkbox"/>	λ	Blue	0	2000	mAU	Cation Exchange	6/13/2012 12:12 PM	Undefined
<input type="checkbox"/>	Conductivity	Red	0	500	mS/cm	Cation Exchange	6/13/2012 12:12 PM	Undefined
<input type="checkbox"/>	%B	Black	0	110	%	Cation Exchange	6/13/2012 12:12 PM	Undefined

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
- λ - Run1
 - λ - Run3
 - %B - Run1
- Show Fractions
 Show Events
 Show Errors

Peak Integration

Auto Integration Manual Integration

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
NGC system: Configuraciones

The background features a light blue gradient on the left and a white-to-yellow gradient on the right. A grid of thin, light blue lines is visible. On the right side, there are several thick, overlapping vertical and horizontal bars in blue, orange, and green. A dark brown square is positioned at the intersection of a blue horizontal bar and an orange vertical bar.

Información para pedidos

<https://www.bio-rad.com/prd/en/US/LSR/PDP/MFCVAQ4VY/NGC-10-ml-Medium-Pressure-Chromatography-Systems>

Home > Life Science Research > Products > Chromatography > Chromatography Systems, Components, and Accessories > NGC™ Medium-Pressure Chromatography Systems > NGC™ 10 ml Medium-Pressure Chromatography Systems

NGC™ 10 ml Medium-Pressure Chromatography Systems

Overview

Configurator

System Tour

NGC 10 ml System Modules

Description Ordering Accessories Fittings Modules Documents

NGC Quest™ 10 Chromatography System 788-0001

Includes 10 ml/min pumps, single-wavelength (UV) and conductivity detection, and automated sample injection, for the purification of proteins

List Price: [Log in](#)

Your Price: [Log in](#)

Quantity:

[Add to Cart](#)

[Add to Hot List](#)

[Add to Quote](#)

Find Products & Information

[Search](#)

Get Help

Call us at 1-800-4-BIORAD (1-800-424-6723)

- » [Email a quick question](#)
- » [Find local contacts](#)
- » [Get product support](#)
- » [Request a quote](#)

Related Products

- » [Design Your NGC™ System](#)
- » [NGC™ Detector Modules and Components](#)
- » [NGC™ Valve Modules and Components](#)
- » [NGC™ Air Sensor Module and Components](#)
- » [NGC™ Pump and Mixer Modules and Components](#)
- » [ChromLab™ Software](#)

News & Events

News

- » [Bio-Rad Receives Frost & Sullivan 2012 North American Laboratory Researchers' Choice: Future Market Leader of Digital PCR Technology Award](#)
- » [Bio-Rad Introduces New Bio-Plex Pro™ RBM Kidney Toxicity Panels](#)

[See all »](#)

Información para pedidos

BUILD YOUR OWN NGC.

The NGC™ chromatography system introduces a new level of modularity and usability for protein purification.

Use our interactive NGC configurator to design your ideal chromatography system.

USE THE PAGE ARROWS TO GET STARTED

1. SELECT A MODEL

Choose an NGC model that best fits your research.

2. DESIGN YOUR SYSTEM

Customize your NGC system with additional modules tailored for your applications.

3. VIEW YOUR NGC

View your personalized NGC system and learn about its capabilities in the Design Summary.

Información para pedidos

1 SELECT A MODEL - 2 - 3

NGC QUEST SYSTEM

NGC QUEST PLUS SYSTEM

NGC SCOUT SYSTEM

NGC SCOUT PLUS SYSTEM

NGC DISCOVER SYSTEM

NGC DISCOVER PRO SYSTEM

THE NGC QUEST™ SYSTEM

Workhorse system for the purification of proteins, peptides, and nucleic acids

- Great for any chromatography lab environment
- Available expansion slots for additional modules
- NGC™ Quest Plus system upgrades to a multi-wavelength detector for increased detection flexibility

THE NGC SCOUT™ SYSTEM

All-purpose purification system with buffer scouting

- Great for more flexible gradient scouting options
- Available expansion tiers for additional modules
- NGC™ Scout Plus system upgrades to a multi-wavelength detector for increased detection flexibility

THE NGC DISCOVER™ SYSTEM

Versatile automated purification system for advanced method development

- Great for large sample numbers and volumes that require more automation and throughput
- Available expansion slots for additional modules
- NGC™ Discover Pro system includes an additional sample inlet valve and an

Información para pedidos

NGC Medium-Pressure Chromatography Systems

NGC Quest Chromatography Systems

For the all-purpose purification of biomolecules

- 788-0001 NGC Quest 10 System
- 788-0003 NGC Quest 10 Plus System
- 788-0002 NGC Quest 100 System
- 788-0004 NGC Quest 100 Plus System

NGC Scout Chromatography Systems

For rapid scouting of biomolecules

- 788-0005 NGC Scout 10 System
- 788-0006 NGC Scout 100 System

For rapid scouting of proteins, peptides, and nucleic acids

- 788-0007 NGC Scout 10 Plus System
- 788-0008 NGC Scout 100 Plus System

NGC Discover Chromatography Systems

For method development

- 788-0001 NGC Discover 10 System
- 788-0003 NGC Discover 100 System

Información para pedidos

Valves

788-4007 **NGC Sample Inject Valve**, pkg of 1, kit includes necessary tubing, fittings, and sample injection port, for manual sample application of small-volume samples via sample loops or large-volume samples using a sample pump

788-4010 **NGC Buffer Blending Valve Module**, pkg of 1, kit includes necessary tubing and fittings, for inline buffer preparation and generating pH gradients for quick pH scouting

788-4006 **NGC Inlet Valve Module**, pkg of 1, kit includes necessary tubing and fittings, for automated switching between multiple buffers and samples during method development

788-4012 **NGC Column Switching Valve Module (10 ml)**, kit includes the necessary tubing and fittings to accommodate the most common column types, holds 5 columns, for use with F10 systems and multiple columns for quick column scouting and reverse flow

788-4026 **NGC Column Switching Valve Module (100 ml)**, kit includes the necessary tubing and fittings to accommodate the most common column types, holds 5 columns, for use with F100 systems and multiple columns for quick column scouting and reverse flow

Air Sensor

788-4017 **NGC Air Sensor Module**, pkg of 1, kit includes 2 large-bore air sensors to protect against air entering pumps and columns, supports up to 4 large- and small-bore air

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